



Digitcool^α: a new performant and ergonomic controlled-rate freezer



Perrine Nogues¹, Agnès Camus¹, Amélia Goanec¹, Charlène Rouillon¹, Lucie Gavin-Plagne¹, Guy Delhomme¹, Eric Schmitt¹

¹IMV Technologies, Z.I. n°1 Est, 61300 Saint-Ouen-Sur-Iton, France
www.imv-technologies.com

+33 (0)672434358 perrine.nogues@imv-technologies.com




CONTEXT

Cryopreservation is an integral part of Assisted Reproduction Technologies in dairy farming. Cryopreservation must be optimized in order to preserve the spermatozoa integrity. IMV Technologies developed a new cryogenic freezer: **Digitcool α**.


OBJECTIVE

This study aimed to compare the quality of bovine semen cryo-frozen in Digitcool α and Digitcool 5300.

MATERIALS & METHODS



RESULTS



CONCLUSION

3

Semen quality was significantly better preserved in Digitcool α. This improvement was observed for all freezing cycles, showing the lowest coefficient of variation.

Digitcool α showed significantly less variability in semen quality between positions. This could be explained by a better thermal exchange homogeneity in the chamber of Digitcool α. Further analyses such as thermal recording in the whole chamber should be performed in order to confirm the better performance of Digitcool α.